



#### **Standard Test Methods**

TAILE S

ASTM G76 - Standard Test Method for Conducting Erosion Tests by Solid Particle Impingement Using Gas Jets **The Falex Air Jet Erosion Test System** controls a repeated impact erosion approach involving a small nozzle delivering a stream of gas containing abrasive particles which impact the surface of a test specimen. The amount of weight loss per unit of time represents the erosion or wear rate, providing laboratory scale erosion measurements under a range of conditions. The Test System may be used to rank the erosion resistance of materials, coatings, platings, epoxies, weldments, sprays, paints, and similar products under the controlled conditions of the test.

The test system uses dry air (100 psi max) and a nozzle tube to control particle velocity and an adjustable media feed nozzle to control the particle feed rate of the abrasive impacting the test sample. An adjustable test sample holder controls the angle of incidence from 15° to 90°.

The test system accommodates standard or user selected abrasive media. The distance from the nozzle to the flat specimen surface is adjustable and is controlled with a micrometer. The test system includes a Start/Stop time controller and a mechanical jet flow shutter, and an exhaust vent (exhaust system not included).

Falex Corporation has specialized in the manufacture of friction, wear and abrasive test equipment for over 75 years.



## **Specifications & Features**

Media Feed Digital display of flow wheel RMP to control the

destiny of sand media. Four pound media feed

chamber.

**Nozzle Pressure** Digital display of pressure at nozzle.

Flow Rate Includes air flow indicator and regulator (Range

0 - 30 LPM).

**Test Time** Digital timer/controller for automatic test

duration.

**Angular Adjustments** Precision scaled adjustment for setting the angle

of incidence of sand against specimen (Range

15-90 degrees).

Nozzle Distance Micrometer adjustment for accurately setting
Adjustment nozzle to specimen distance (Range 0.075 to

1.000 inch).

Vacuum Exhaust Minimum of one inch of mercury column

required (source not provided).

**Pneumatic** 0-100 psig required (minimum pipe inside

diameter 3/8") (source not provided).

**Electrical** 220 volts AC single phase, with a

15 amperage requirement.

**Wear** Measured by the amount of weight loss per unit

of time. This represents the erosion or wear rate.

**ASTM Standard** G76 - Standard Test Method for Conducting

Erosion Tests by Solid Particle Impingement

using Gas Jets

# **Ordering Information**

## Part Number Description

007-001-001	Falex Air Jet Erosion Test
	Machine Complete

### Reference Material (ASTM G76) Test Specimens

007-501-001	Test Panel, 1020 Steel, HRB 68-72, 30 mm x 10 mm x 2 mm
007-502-002	Test Panel, 1020 Steel, HRB 68-72, 2" x 2" x 1/16"
007-599-001	Aluminum Oxide Test Abrasive, Grade 240 grit, 5 lb.

Falex Corporation follows a policy of continuous product improvement. Specifications are subject to change without notice.



### Distributed by:

